

Luer slip IV syringes




**IVS01** 1ml Luer Slip IV Syringe (100)



**IVS01INS** 1ml Luer Slip Insulin Syringe (100)



**IVS03** 3ml Luer Slip IV Syringe (100)



**IVS05** 5ml Luer Slip IV Syringe (100)



**IVS10** 10ml Luer Slip IV Syringe (100)




**IVS10E** 10ml Eccentric Luer Slip IV Syringe (100)



**IVS20** 20ml Luer Slip IV Syringe (50)




**IVS30** 30ml Luer Slip IV Syringe (30)



**IVS60** 60ml Luer Slip IV Syringe (30)


Luer lock IV syringes




**IVL01** 1ml Luer Lock IV Syringe (100)




**IVL03** 3ml Luer Lock IV Syringe (100)



**IVL05** 5ml Luer Lock IV Syringe (100)




**IVL10** 10ml Luer Lock IV Syringe (100)



**IVL20** 20ml Luer Lock IV Syringe (50)




**IVL30** 30ml Luer Lock IV Syringe (30)




**IVL60** 60ml Luer Lock IV Syringe (30)

Luer lock insulin syringe




**IVL60INS** 60ml Luer Lock Insulin Syringe (30)

disposable tourniquets



**DT02** Adult Disposable Tourniquet (100)



**DT03** Paediatric Disposable Tourniquet (100)

blunt fill needles



**BN01** Blunt Fill Needle (100)



**BN02** Blunt Fill Needle with Filter (100)



# IV syringes



Medicina  
IV syringes

Cost saving solutions  
without compromising  
safety and efficacy



Units 1-4, Rivington View Business Park,  
Station Road, Blackrod, Bolton BL6 5BN, UK  
Telephone: +44(0)1204 695050  
[www.medicina.co.uk](http://www.medicina.co.uk)



[www.facebook.com/MedicinaUK](http://www.facebook.com/MedicinaUK)



@medicina\_ltd

leaders in enteral delivery





# Medicina Luer slip IV syringes

Medicina has manufactured a complete range of Luer slip IV syringes that can be used with any standard, filtered or safety needle. The overall aim is to reduce the cost to the healthcare sector whilst improving patient safety.

The Luer slip syringe is available in a range of barrel sizes:

- 1ml
- 3ml
- 5ml
- 10ml
- 20ml
- 30ml
- 60ml

Sizes 10ml through to and including 60ml have an eccentric tip. This gives the end user a clear view of the entry site and allows the site to be entered at the correct angle.

The flat plunger fits flush against the flat shoulder of the barrel, eliminating dead space to ensure accurate measurements are delivered to the patient.



Medicina have printed the wording "IV Syringe" onto the barrel of each size of syringe, making the use of this product clear, therefore avoiding any confusion and improving patient safety.

# Medicina Luer lock IV syringes

Medicina has developed a range of Luer lock IV syringes for use with IV extension sets, standard or safety needles and IV cannula. They are designed to be used on the most common types of syringe drivers.

The clear print on the barrel of the Medicina Luer lock IV syringe ensures accuracy when drawing up fluid. The syringe is clearly labelled as an IV syringe, making its use clear and avoiding any confusion.

The Luer lock syringe is available in a range of barrel sizes:

- 1ml
- 3ml
- 5ml
- 10ml
- 20ml
- 30ml
- 60ml

The Medicina syringe also eliminates dead space as the cone shaped plunger fits flush against the shoulder of the barrel of the syringe. This further ensures that the patient receives the correct amount of fluid.

The Medicina Luer lock syringes have been designed to be used with the most common brand of IV syringe drivers. Medicina also offer a 60ml Luer lock insulin syringe. The red print on the barrel clearly displays that the syringe is "FOR INSULIN ONLY". This is aimed to improve patient safety.

Improving safety whilst cutting costs.

